



Hubbadewanne AVERO Phoenix
Lifting tub AVERO Phoenix 170/190

Art. Nr. 910130000 (Typ 170)
Art. Nr. 910131000 (Typ 190)

Bauseitige Anschlüsse - Sanitär - Elektro

1. Warmwasseranschluß ¾ Zoll, Innengewinde Plattenbündig
2. Kaltwasseranschluß ¾ Zoll, Innengewinde Plattenbündig
3. Abfluß DN 50, plattenbündig
4. Abzweigdose IP 65, 230V/50Hz, 1,5 mm², Absicherung 10 A.

Bauseitige Schutzmaßnahmen

Vorabspernung Kalt- und Warmwasser, Allpoliger Ausschalter innerhalb des Raumes, Fi-Schutzschalter - Nenn-Fehlerstrom max. 30mA, Elektroinstallation nach den zur Zeit gültigen DIN/VDE Normen. Potentialausgleich nach DIN VDE 0107.

Pre-installed connections - Sanitary - Electric

1. Warm water connection ¾", Inside thread, tile-even
2. Cold-water connection ¾", Inside thread, tile-even
3. Drainage DN 50, 50 mm above tiles/floor
4. Junction-box Ip65, 230V/50Hz, 1,5 mm², Protection 10 Ampere.

Pre-installed security measure

Pre-closing valves cold-/warm water.
Electric connections must be done according to the rules and regulations for electric used devices in wet-rooms for each country.